

NEUROCRYOSTIMULATION: -78°C° IN YOUR HANDS

CRYOSCREEN *



RELIEVE PAIN QUICKLY WITH NEUROCRYOSTIMULATION (NCS).

This technique makes it possible to drop the temperature from 32°C to 4°C in 30 seconds thanks to a cold power of -78°C and 50 bars of pressure (CO₂).

THE THERMAL SHOCK IS IMMEDIATE

It stimulates the skin receptors that send a powerful and very fast message to the brain to get a response by the reflex system.

The results are immediately felt by patients as:

- * A pain relief effect
- * Management of the inflammation crisis
- * Muscular relaxation
- * A vasomotor effect that allows rapid draining of hematomas and edemas



Some examples of pathologies treated:

- 1 Joint problems with inflammatory flare-ups
- 2 Epicondylitis
- 3 Tendinopathies
- 4 The aftermath of muscular accidents
- 5 Post-traumatic hematomas and edemas

MASTERY OF -78°C IN HIGHLY CONCENTRATED TECHNOLOGY

Using an embedded thermal camera and powerful screen, you measure the temperature from a distance and visualise the heat zones to be treated.

A real tool to rapidly detect, target and take action for better patient comfort. Its ergonomics, its size and its handling allow it to be used in any situation and everywhere.



Performance is made possible thanks to precise measurement, quality imaging and the ability to take photos of the area before, during and after treatment.

Control of the thermal shock is assured as well as visualisation of the results.

NeuroCryoStimulation is now part of the therapeutic arsenal that supports the practitioner's action and the patient every day for treatment in acute or chronic phases.

www.cryonic-medical.com

Cryonic Medical SAS 6 rue Berthelot, 25000 Besançon - France
+33 3 81 54 35 40 / cryonic@cryonic-medical.com / www.cryonic-medical.com



CRYONIC
NEUROCRYOSTIMULATION
COLD TREATMENT TECHNIQUE

NEUROCRYOSTIMULATION: -78°C° IN YOUR HANDS

CRYOSCREEN *



RELIEVE PAIN QUICKLY WITH NEUROCRYOSTIMULATION (NCS).



This technique makes it possible to drop the temperature from 32°C to 4°C in 30 seconds thanks to a cold power of -78°C and 50 bars of pressure (CO₂).

THE THERMAL SHOCK IS IMMEDIATE

It stimulates the skin receptors that send a powerful and very fast message to the brain to get a response by the reflex system.

The results are immediately felt by animal as:

- * A pain relief effect
- * Management of the inflammation crisis
- * Muscular relaxation
- * A vasomotor effect that allows rapid draining of hematomas and edemas



Some examples of pathologies treated:

- 1 Tendinopathy (lameness)
- 2 Fetlock sprains
- 3 Ligament inflammations
- 4 Haematomas and oedema
- 5 Muscular spasms (back pain)

MASTERY OF -78°C IN HIGHLY CONCENTRATED TECHNOLOGY



Using an embedded thermal camera and powerful screen, you measure the temperature from a distance and visualise the heat zones to be treated.

A real tool to rapidly detect, target and take action for better animal comfort. Its ergonomics, its size and its handling allow it to be used in any situation and everywhere.

Performance is made possible thanks to precise measurement, quality imaging and the ability to take photos of the area before, during and after treatment.

Control of the thermal shock is assured as well as visualisation of the results.

NeuroCryoStimulation is now part of the therapeutic arsenal that supports the daily care for acute or chronic pathologies.



www.cryonic-medical.com

Cryonic Medical SAS 6 rue Berthelot, 25000 Besançon - France
+33 3 81 54 35 40 / cryonic@cryonic-medical.com / www.cryonic-medical.com

